

# Zapatos para Diabéticos de NursingCare

Medidas del pie en mm para diferentes tipos de zapatos.

	X	Y	Z	
35	227	96	234	230
36	233	97	239	235
37	239	98	243	239
38	246	100	248	244
39	253	101	252	248
40	260	103	257	253
41	267	105	262	258
42	273	106	266	262

	X	Y	Z	
39	252	104	252	248
40	259	105	257	253
41	265	106	261	257
42	272	107	265	262
43	278	108	270	266
44	285	109	275	271
45	291	111	279	275
46	298	112	284	280

	X	Y	Z	
39	252	105	270	266
40	259	106	275	271
41	265	107	279	275
42	272	108	284	280
43	278	110	288	284
44	285	111	293	289
45	291	113	297	293
46	298	114	302	298

	X	Y	Z	
39	252	108	266	262
40	259	109	271	267
41	265	110	275	271
42	272	111	280	276
43	278	113	284	280
44	285	114	289	285
45	291	115	293	289
46	298	116	298	294

	X	Y	Z	
39	252	105	256	252
40	259	106	261	257
41	265	107	265	261
42	272	107	270	266
43	278	109	274	270
44	285	110	279	275
45	291	111	283	279
46	298	112	288	284

	X	Y	Z	
39	252	101	246	242
40	259	102	251	247
41	265	103	255	251
42	272	103	260	256
43	278	104	264	260
44	285	105	269	265
45	291	106	273	269
46	298	107	278	274



Talla



Longitud del pie



Ancho del pie

**Z**

Colocación de pies



Sinplantilla de media corcho.



Con plantilla de media corcho.

## Forma de pie

**A**

Extra ancho

**B**

Ancho

**C**

Normal

